



EMAR Beijing Emarbio Science & Technology Co.,LTD

anti-Bax Rabbit Polyclonal antibody

Catalog No.	EM51397-01	EM51397-02
Size	50 μ l	100 μ l
Species reactivity:	Human Mouse Rat	Molecular Weight: 21KD

All Names:Apoptosis regulator BAX; BAX; Bcl-2-like protein 4; BCL2-associated X protein; Bcl2-L-4; BCL2L4

Immunogen:A synthesized peptide derived from human Bax

Purification:Affinity-chromatography


Application: WB 1:500-1:3000 IHC 1:50-1:200 IF 1:200

Optimal dilutions/concentrations should be determined by the end user.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis on HepG2 cell lysate using Bax Antibody

EMAR

EMAR