



Anti-GRB2 Rabbit Polyclonal Antibody

Catalog No.	EM55472-01	EM55472-02
Size	50 μ l	100 μ l
Species reactivity:	Human Mouse Rat	Molecular Weight: 32KD

Background:

All Names: Growth factor receptor-bound protein 2; Adapter protein GRB2; SH2/SH3 adapter GRB2;

Immunogen: A synthesized peptide

Specificity: GRB2 Antibody detects endogenous levels of total GRB2


Application: WB 1:1000-1:3000

Optimal dilutions/concentrations should be determined by the end user.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of extracts from Jurkat cells, using GRB2 antibody.