



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-CBP Tag Mouse Monoclonal Antibody(Mix)

Catalog No.	EM33208-01	EM33208-02
Size	50 μ l	100 μ l

Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.

Molecular Weight: N/A

Application : WB 1:5,000-10,000 IP 1:200 Optimal dilutions should be determined by the end user.

Specificity:The CBP tag antibody can detects CBP recombinant protein.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of His-MBP-V5 recombinant protein with V5 mAb diluted at 1:10,000.

EMAR

EMAR