



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-TAP Tag Mouse Monoclonal Antibody(Mix)

Catalog No.	EM33209-01	EM33209-02
Size	50μl	100μl

Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.

Molecular Weight: N/A

Application : WB 1:5,000-10,000 IP 1:200 Optimal dilutions should be determined by the end user.

Specificity:The TAP tag antibody can detects TAP recombinant protein.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of TAP recombinant protein with TAP mAb diluted at 1:5,000.

EMAR

EMAR