



EMAR Beijing Emarbio Science & Technology Co.,LTD

Anti-EFHD1 Rabbit Polyclonal Antibody

Catalog No.	EM42024-01	EM42024-02
Size	50μl	100μl
Species	Human	Molecular Weight: 27KD

Background: EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation.

Specificity: The EFHD1 Rabbit Polyclonal antibody detects endogenous EFHD1 proteins.


Application: WB: 1:1,000

Optimal dilutions/concentrations should be determined by the end user.

Form: Liquid ,1mg/ml

Storage buffer: PBS, pH 7.4, containing 0.02% **sodium azide** as Preservative and 50% Glycerol.

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of 293T cell lysate with EFHD1 rabbit polyclonal antibody diluted at 1:2,000.

EMAR

EMAR