



anti-iNOS Rabbit Polyclonal antibody

Catalog No.	EM52262-01	EM52262-02
Size	50µl	100µl
Species reactivity:	Human	Molecular Weight:

All Names: NITRIC OXIDE SYNTHASE, INDUCIBLE, NITRIC OXIDE SYNTHASE, MACROPHAGE, NOS2A, INDUCIBLE, HEPATOCYTE, nitric oxide synthase 2A (inducible, hepatocytes)

Immunogen: The antiserum was produced against synthesized peptide derived from human iNOS.


Purification: The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Application: WB 1:500-1:3000 IHC 1:50-1:200 IF 1:100-1:500 ELISA 1:20000
Optimal dilutions/concentrations should be determined by the end user.


Form: Liquid ,1mg/ml

Storage buffer: Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. Store at -20 ° C. Stable for 12 months from date of receipt

Storage: Store at -20°C. Do not aliquot the antibody.



Western blot analysis of extracts from NIH/3T3 cells, using iNOS antibody



Immunohistochemical analysis of paraffin-embedded human lung carcinoma tissue, using iNOS antibody

